

Cosmopore®

Soft, absorbent latex-free adhesive dressing for superior wound management.



Cosmopore®

Latex-Free Adhesive Wound Dressing

Cosmopore® provides a sterile wound contact layer that is air and water permeable, 100% cotton and non-adherent. This pad performs well to absorb exudate for a cleaner wound bed. The non-woven adhesive cover of the Cosmopore pad has rounded edges to prevent the dressing from becoming detached even after being worn for prolonged periods of time. The cover is coated with a non-sensitizing polyacrylate adhesive that is well tolerated by the skin.

Description

- ▶ Latex-free, sterile, individually sealed
- ▶ 100% cotton non-adherent island pad absorbs exudate for a cleaner wound bed
- ▶ Non-woven adhesive cover has rounded edges to prevent detachment
- ▶ Wide adhesive border ensures excellent adhesion
- ▶ Convenient, easy-to-apply non-adherent adhesive dressings for Athletic Trainers

Applications

- For use as a post-operative dressing and for minor wounds
- Square dressings can be used as cover dressings for alginates such as Sorbalgon® and gel sheets such as AquaClear®



Cosmopore® adapts well to body contours and is highly absorbent.

Cosmopore®		Dressings	Dressings
Order No.	Size	Per Box	Per Case
48000000	2" x 3 1/2"	50	400
48100000	3 1/8" x 4"	25	200
48200000	3 1/8" x 6"	25	200
48300000	4" x 8"	25	200
48400000	4" x 10"	25	200
48500000	4" x 13 3/4"	25	200
49500000	4" x 4"	25	200
49510000	6" x 6"	25	200